

AJAX "VIOLET RAY" GENERATOR
 (High Frequency)

OPERATING INSTRUCTIONS

(1) Select the Electrode it is intended to use and insert gently into the handle. Do not FORCE it; the seating of square end is easily found.

(2) Switch off "Mains" Current and connect the two pin plug or adaptor to the electric supply.

(3) Push the sliding switch on apparatus as far as it will go towards "OFF".

(4) Place the electrode firmly on the part to be treated and then switch on the current at the Mains.

To increase the strength push the sliding switch forward towards "Mild", "Medium" and "Strong" or for Lumbago and deep seated ailments the treatment can be given through an article of clothing instead of direct on to the skin.

(5) When the treatment is finished bring sliding switch back to "OFF" before losing the contact of the electrode with the part treated and then switch off current at Mains.

The above can be followed quite easily by referring to the illustrations in the booklet.

(6) The lubricant supplied with the apparatus should be applied very sparingly to the part to be treated or alternatively a small drop can be placed on the electrode.

A fresh supply (full size bottle) can be obtained from us, price 2/6d.

(7) In treating the surface of the body or scalp always keep the Electrode moving slowly.

(8) Before treating the scalp, it should be thoroughly cleansed of all grease, lotions, hair tonics and other preparations.

(9) Keep the electrodes clean.

(10) Never insert or remove an electrode with the current on.

(11) Never use the apparatus when in the bath or when standing on a wet floor.

(12) Do not leave current switched on when not in use.